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**DEI® Experience Evolution DENTINA PASTA** (1 siringa da 5 g):  
A1/B1 (REF. EEP1), A2/B2 (REF. EEP2), A3 (REF. EEP3),  
A3.5 (REF. EEP4), A4 (REF. EEP5), B3/B4 (REF. EEP6),  
C1/D2/D3/D4 (REF. EEP7), C2/C3/C4 (REF. EEP8)  
**DEI® Experience Evolution DENTINA FLOW** (1 siringa da 3 g):  
A1/B1 (REF. EEF1), A2/B2 (REF. EEF6), A3 (REF. EEF2),  
A3.5 (REF. EEF7), A4 (REF. EEF3), B3/B4 (REF. EEF4),  
C1/D2/D3/D4 (REF. EEF8), C2/C3/C4 (REF. EEF5)  
**DEI® Experience Evolution SMALTO** (1 siringa da 5 g):  
LIGHT (REF. EXSL), MEDIUM (REF. EXSM), DARK (REF. EXSD)  
**DEI® Experience Evolution INCISALE** (1 siringa da 5 g):  
BLUE (REF. EEIB), CLEAR (REF. EEIC), LIGHT (REF. EEIL)  
**DEI® Experience Evolution STAIN** (1 siringa da 1 g):  
ORANGE (REF. EESO), RED (REF. EESR), DARK BROWN (REF. EESDB)

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THE NEW COMPOSITE SYSTEM FOR  
UNIVERSAL RESTORATIONS  
AND PREPARATIONS OF BRIDGES,  
CROWNS, INLAYS, FACETS,  
AND CHARACTERIZATION  
OF MOBILE PROSTHESES

Product for dental use  
Keep away from children  
Date of the specifications: 04/2021

**AVAILABLE PACKAGES**

**DEI® Experience Evolution DENTINA PASTA** (One 5 g syringe):  
A1/B1 (REF. EEP1), A2/B2 (REF. EEP2), A3 (REF. EEP3),  
A3.5 (REF. EEP4), A4 (REF. EEP5), B3/B4 (REF. EEP6),  
C1/D2/D3/D4 (REF. EEP7), C2/C3/C4 (REF. EEP8)  
**DEI® Experience Evolution DENTINA FLOW** (One 3 g syringe):  
A1/B1 (REF. EEF1), A2/B2 (REF. EEF6), A3 (REF. EEF2),  
A3.5 (REF. EEF7), A4 (REF. EEF3), B3/B4 (REF. EEF4),  
C1/D2/D3/D4 (REF. EEF8), C2/C3/C4 (REF. EEF5)  
**DEI® Experience Evolution SMALTO** (One 5 g syringe):  
LIGHT (REF. EXSL), MEDIUM (REF. EXSM), DARK (REF. EXSD)  
**DEI® Experience Evolution INCISALE** (One 5 g syringe):  
BLUE (REF. EEIB), CLEAR (REF. EEIC), LIGHT (REF. EEIL)  
**DEI® Experience Evolution STAIN** (One 1 g syringe):  
ORANGE (REF. EESO), RED (REF. EESR), DARK BROWN (REF. EESDB)

**CHARACTERISTICS**

It was very difficult to achieve the evolution of a material like DEI® experience which has become a reference for all composites. We worked on three aspects:  
1. Simplification of shades with a chameleon effect, 8 dentine paste, 8 dentine flow and 3 enamels are enough to obtain all the shades of the Vita® scale.  
2. Improvement of technical data, strength, modulus of elasticity, water absorption, abrasion.  
3. We have added the Dentina Flow for large first and second class cavities in the studio and for molding in the laboratory which we will use as the core of the elements.  
The consistency of Dentina Flow will significantly increase adhesion with the tooth and with the substrate thanks to greater contact / wettability, and will avoid, in the molding technique, joint elevation after pressing in the flask.  
Using Dentine Flow inside, with dentin, enamel and incisal paste, we will obtain the consistency of a natural tooth, elastic inside, very hard, with low abrasion of the enamel on the outside and greater shock absorption capacity.  
As regards the new M.C.M.® technique, we have created an explanatory video tutorial (www.deiitalia.it).  
To achieve the shades, consult the outline of use of the mass for the correct use of DEI® Experience Evolution.  
Water absorption: 40 µg/mm³

**COMPOSITION**

**Composition Dentina Pasta**  
- Anorganic filler (74% by weight and 66% by volume):  
Glass filler (average particle size: 0.7 µm)

Pyrogenic silica (average particle size: 0.04 µm)  
- Monomers (25% by weight):  
Urethane dimethacrylate, 1,4-butanediol dimethacrylate, BisGMA  
- Additives (1 % by weight): initiators, stabilizers, pigments  
**Composition Dentina Flow**  
- Anorganic filler (57% by weight, 51% by volume)  
Glass filler (average particle size: 0.7 µm)  
Pyrogenic silica (average particle size: 0.04 µm)  
- Monomers (43 % by weight):  
Urethane dimethacrylate, 1,4-butanediol dimethacrylate  
- Additives (1% by weight): initiators, stabilizers, pigments  
**Composition Smalto**  
- Anorganic filler (75% by weight, 53% by volume) (0,005 - 3,0 µm)  
- Glass powder, diurethane dimethacrylate, silicon dioxide, Bis-GMA, tetramethylene dimethacrylate.  
**Composition Incisale**  
- Anorganic filler (75% by weight, 52% in volume):  
Glass filler (average particle size: 0,7µm)  
Pyrogenic silica (average particle size: 0,04 µm)  
- Monomers: Diurethane dimethacrylate, 1,4-butanediol dimethacrylate, Bis-GMA  
**Composition Stain**  
- Anorganic filler (50% by weight 36% by volume):  
Glass filler (average particle size: 0.7 µm)  
Pyrogenic silica (average particle size: 0.01 µm)  
- Monomers (49% by weight):  
Urethane dimethacrylate, Tetramethylene dimethacrylate, Bis-GMA  
- Additives (1% by weight): catalysts, stabilizers, pigments

**WAY OF USE (STUDIO)**

Proceed as usual preparing the cavity, protecting the pulp, etching and applying the adhesive system. Use of the adhesive system DEI® New Extra Bond One (not included in the kit) is recommended. Apply the composite material stratifying the various masses. Every layer must be cured for 40 seconds for depths not thicker than about 1.5 mm or cured for 10 sec. with DEI® experience LED Studio. Photo-polymerize with curing lamps with light performance of at least 500 mW/cm2. If it is impossible to put the light source close enough to the material, it is necessary to lengthen the exposure time.

The observance of photo-polymerization times and depths are fundamental for achieving complete conversion of the materials, and ensuring resistance and duration in the time. The composites harden also very few seconds, but are not polymerized. In this phase saving time could cause failure of the restoration.

**ADHESION ON CURED COMPOSITE**

To get perfect adhesion between already cured and finished composite and composite not yet cured (for example in case of repairs): roughen a little the already cured composite, clean it with ethylic alcohol and clean brush; apply on the so prepared surface a thin layer of DEI® experience Seal Coat Fast; don't cure and add new composite. Light cure for 60 seconds with halogen or led lamp for dentist surgeries (or for 8 minutes with DEI® experience LED in laboratory).

**FINISHING TECHNIQUE (STUDIO)**

Before light curing the last layer of composite, apply DEI® Surface Hardener and light cure for 40 seconds, wash and dry. The hardening of the composite shall be complete. The use of that product avoids the need of eliminating the surface inhibition layer with burs or other items, consequently with shorter working times. After modeling, run a micro-brush with ethylic alcohol on the surface. With disposable brush spread a thin layer of DEI® Seal Coat Fast Led on the surface, blow with clean and oil/water free air to improve the distribution of the product; run a clean brush or micro-brush for elimination of surpluses, let work 10 seconds and light cure for 90 seconds. Take out the surface layer of DEI® Seal Coat Fast Led using DEI® Polish Paste and goat hair brush, and polish with cotton small brush. DEI® Seal Coat Fast Led remains sealing the micro-porosity. At the end of polishing an extremely hard, resistant and shiny surface shall be obtained. In addition, the present technique allows considerable diminution of growth of tartar, bacterial plaque and of colorations caused by beverages such as coffee.

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**FORMAS DE SUMINISTRO**

**DEI® Experience Evolution DENTINA PASTA** (1 jeringa de 5 g):  
A1/B1 (REF. EEP1), A2/B2 (REF. EEP2), A3 (REF. EEP3),  
A3.5 (REF. EEP4), A4 (REF. EEP5), B3/B4 (REF. EEP6),  
C1/D2/D3/D4 (REF. EEP7), C2/C3/C4 (REF. EEP8)  
**DEI® Experience Evolution DENTINA FLOW** (1 jeringa de 3 g):  
A1/B1 (REF. EEF1), A2/B2 (REF. EEF6), A3 (REF. EEF2),  
A3.5 (REF. EEF7), A4 (REF. EEF3), B3/B4 (REF. EEF4),  
C1/D2/D3/D4 (REF. EEF8), C2/C3/C4 (REF. EEF5)  
**DEI® Experience Evolution SMALTO** (1 jeringa de 5 g):  
LIGHT (REF. EXSL), MEDIUM (REF. EXSM), DARK (REF. EXSD)  
**DEI® Experience Evolution INCISALE** (1 jeringa de 5 g):  
BLUE (REF. EEIB), CLEAR (REF. EEIC), LIGHT (REF. EEIL)  
**DEI® Experience Evolution STAIN** (1 jeringa de 1 g):  
ORANGE (REF. EESO), RED (REF. EESR), DARK BROWN (REF. EESDB)

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**WAY OF USE (LAB)**

DEI® experience is the suitable composite for the application with pressing technique, known as M.C.M.® (Monolithic Composite Method). We suggest to dental technicians to join us in free work shops held by the Odt. Paolo Pagliari (creator of the technique) at our headquarters, or watch the video tutorials on www.deiitalia.it.

**GUARANTEE**

Our technical instructions, regardless of whether they are provided verbally, in writing or during practical demonstrations, are based on our own experience and should only be considered guidelines. As our products are subject to continued development, we reserve the right to modify them.

**STORAGE**

To ensure the product stability, store the material between 10°C and 25°C. If, during high temperature periods outside, the product is kept in a refrigerator, it must be removed well before use to resume its original contouring properties. Close the syringe tightly after use and protect it against direct light. Turn back spindle about one turn, in order to keep the material from leaking out. If the material is not cured fully, it may discolor, its full mechanical properties will not be achieved and irritation of the pulp cannot be ruled out. This product was developed specifically for the range of indications described. It must be used as described in the instructions. The manufacturer cannot be held responsible for damage due to incorrect handling or processing.

**SHELF-LIFE**

The maximum shelf-life is printed on the label of each syringe. Do not use after the expiry date.

**SIDE-EFFECTS**

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects - even if it is doubtful that the side-effect has been caused by our product - please kindly contact us. To prevent possible reactions of the pulp in cavities where the dentin is exposed, the pulp must be protected adequately (e. g. calcium hydroxide preparation).

**CONTRA-INDICATIONS / INTERACTIONS**

Contains tetramethylene dimethacrylate, diphenyl (2,4,6-trimethylbenzoyl) phosphine oxide

**Warning.** May cause an allergic skin reaction. Avoid breathing vapours / spray. Wear protective gloves / protective clothing / eye protection / face protection. If skin irritation or rash occurs: Get medical advice / attention.

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product.

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**INTERACTIONS WITH OTHER SUBSTANCES**

As phenolic substances (such as eugenol) inhibit polymerization, do not use cavity liners (such as zinc-oxide eugenol cements) containing such substances.

\* Vita is a registered trademark of Vita® Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.

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\* This is a guideline. Each color includes a value: the higher the value, the lighter it is; the lower the value, the darker it is. If, for example, we need to make an A3.5 with a lighter value, we will use the Smalto Light; if darker, we will use the Smalto Dark. Another component of the color is the chroma, if we are to make an element with more chroma, we must have a higher dentine space. The last component of the color is the tint (Vita® shade guide). We recommend, when a color is detected, to take a photograph of the filed tooth (dentin color), one of an adjacent tooth (for the Vita® shade guide) and finally one in black and white (to determine the value, which is the gray of the shade). In all three photos, near the tooth, there must be the Vita® shade sample.

TECHNICAL DATA	DENTINA PASTA	DENTINA FLOW	SMALTO	INCISALE
Vickers hardness (MPa)	400	370	539	560
Flexural strength (MPa)	95	140	120	140
Modulus of elasticity (MPa)	6900	6500	9000	10000
Compressive strength (MPa)	360	360	400	450
Depth of polymerization (mm)	2.5	2.0	3.0	3.0
Filler content by weight	74%	57%	75%	50%
Filler content by volume	66%	51%	53%	36%

PROBLEM	CAUSE	REMEDY
Composite does not cure properly	Light output of the light-curing lamp is inadequate  Emitted wavelength range of the lightcuring lamp is inadequate	Check the light output and change the light source if required  Consult the manufacturer of the lightcuring lamp. Recommended wavelength range: 350 - 500 nm
Composite in the syringe is sticky and soft, colorless liquid separates in the syringe	Material has been stored for a longer period at 25°C (77°F)  Material has been kept in a syringe warmer for too much time	Adhere to storage temperature. Store at 10 - 25°C (50 - 77°F)  Never keep a syringe in a syringe warmer for more than one hour per application
Composite appears too hard and firm in the syringe	Material stored at temperatures 10°C (50°F) for a longer period of time  Syringe not properly sealed, composite partially cured	Allow the composite to heat to room temperature before use, use a syringe warmer if necessary.  Always seal the syringe properly with the cap after taking out composite
Inlay/onlay is not properly retained when fitted	Restoration is too opaque to be cemented using only light-curing composite	Use dual-curing luting composite
Composite does not cure completely (dark or opaque shades)	Composite layers applied too thickly for each curing cycle	Adhere to a max. thickness of 2.0 mm per layer
Restoration appears too yellow compared with the shade guide	Inadequate curing of the composite layer	Repeat the exposure cycle several times; min. 40 sec.

Via Torino, 765 - 21020 Mercallo (Va) - Italy  
Phone +39 0331 969270 - www.deiitalia.it

ESQUEMA DE EMPLEO DE LOS COMPONENTES			
COLOR SHADE VITA*	OPACO	DENTINA PASTA / FLOW	ESMALTO
A1	LIGHT	A1 / B1	LIGHT
A2	LIGHT	A2 / B2	LIGHT
A3	LIGHT	A3	LIGHT
A3.5	MEDIUM	A3.5	MEDIUM
A4	MEDIUM	A4	DARK
B1	LIGHT	A1 / B1	LIGHT
B2	MEDIUM	A2 / B2	LIGHT
B3	LIGHT	B3 / B4	MEDIUM
B4	MEDIUM	B3 / B4	DARK
C1	MEDIUM	C1 / D2 / D3 / D4	LIGHT
C2	MEDIUM	C2 / C3 / C4	LIGHT
C3	DARK	C2 / C3 / C4	MEDIUM
C4	DARK	C2 / C3 / C4	DARK
D2	LIGHT	C1 / D2 / D3 / D4	LIGHT
D3	DARK	C1 / D2 / D3 / D4	LIGHT
D4	DARK	C1 / D2 / D3 / D4	MEDIUM

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\* Vita es una marca registrada de Vita Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen